

EnzyChrom™ Aspartate Transaminase Assay Kit (Cat# EASTR-100)

Colorimetric Determination of Aspartate Transaminase activity

DESCRIPTION

Aspartate Transaminase (AST), also known as serum glutamic oxaloacetic transaminase (GOT) or aspartate aminotransferase (ASAT/AAT), facilitates the conversion of aspartate and α -ketoglutarate to oxaloacetate and glutamate. There are two isoenzymes in humans: GOT1 is a cytosolic isoenzyme derived from red blood cells and heart; GOT2 is the mitochondrial isoenzyme found mainly in the liver. AST is elevated in liver and muscle diseases. It is part of diagnostic tests for liver function, myocardial infarction, acute pancreatitis, acute hemolytic anemia, severe burns, acute renal disease and trauma.

Simple, direct and automation-ready procedures for measuring AST activity find wide applications in research and drug discovery. BioAssay Systems' AST activity assay is based on the quantification of oxaloacetate produced by AST. In this assay, oxaloacetate and NADH are converted to malate and NAD by the enzyme malate dehydrogenase. The decrease in NADH absorbance at 340 nm is proportionate to AST activity.

KEY FEATURES

Sensitive. Linear detection range: 2–100 U/L.

Simple and convenient. This simple, convenient assay can be carried out in a microplate or a cuvette and takes only 10 min.

APPLICATIONS

Direct Assays: AST activity in serum, plasma and other biological samples.

Drug Discovery/Pharmacology: effects of drugs on AST activity.

KIT CONTENTS

Assay Buffer: 24 mL Cofactor: 120 μ L
Enzyme Mix: 120 μ L NADH: 500 μ L

Storage conditions. The kit is shipped on ice. Store all reagents at -20°C . Shelf life of six months after receipt.

Precautions: Reagents are for research use only. Normal precautions for laboratory reagents should be exercised while using the reagents. Please refer to Material Safety Data Sheet for detailed information.

ASSAY PROCEDURES

Equilibrate all components to room temperature. Keep thawed enzymes on ice.

Assays can be performed at 37°C or at room temperature. Prior to assay, bring the working reagents, microplate and spectrophotometer to the desired temperature.

Assay is compatible with serum or plasma (heparin, EDTA). Samples should be clear and free of particles or precipitates. Hemolyzed samples should not be used.

Procedure using 96-well plate

1. Samples and controls. Transfer 20 μ L sample to each well. For each assay plate, include two wells with 20 μ L distilled water to be used for the NADH Standard and Blank. Keep plate at the desired temperature (e.g. 37°C).

2. Prepare Working Reagent for Sample and Standard wells, by mixing for each well, 200 μ L Assay Buffer, 1 μ L Cofactor, 1 μ L Enzyme Mix and 4 μ L NADH. Warm to desired temperature (e.g. 37°C).

Prepare Blank Reagent for the Blank well, by mixing 200 μ L Assay Buffer, 1 μ L Cofactor, 1 μ L Enzyme Mix and 4 μ L H_2O . Warm to desired temperature (e.g. 37°C).

3. Add 200 μ L Working Reagent to the Standard and Sample wells, and 200 μ L Blank Reagent to the Blank well.

Immediately tap plate to mix, incubate at the desired temperature and read $\text{OD}_{340\text{nm}}$ at 5 min and at 10 min. Alternatively, record kinetics at 340 nm.

Procedure using cuvettes

1. For each assay, include one Standard and one Blank control.

For each Sample and Standard, prepare Working Reagent by mixing 1000 μ L Assay Buffer, 5 μ L Cofactor, 5 μ L Enzyme Mix and 20 μ L NADH. Transfer 990 μ L Working Reagent to each sample cuvette and standard cuvette. Warm to desired temperature (e.g. 37°C).

To Blank control cuvette, add 960 μ L Assay Buffer, 5 μ L Cofactor, 5 μ L Enzyme Mix and 20 μ L H_2O . Warm to desired temperature (e.g. 37°C).

2. Prewarm sample to the desired temperature. Add 100 μ L Sample to the Sample Cuvette. Transfer 100 μ L H_2O to the Standard cuvette and to Blank Control cuvette, respectively.

3. Mix immediately. Read $\text{OD}_{340\text{nm}}$ at 5 min and 10 min. Alternatively, record kinetics at 340 nm.

CALCULATION

For each Sample, calculate the rate of NADH consumption by subtracting the OD at 10 min from the OD at 5 min (ΔOD_S). Similarly, calculate the rate ($\Delta\text{OD}_{\text{NADH}}$) for the NADH standard ($\text{OD}_{5\text{ min}} - \text{OD}_{10\text{ min}}$).

Determine AST activity using the following equation,

$$\text{AST} = 388 \times \frac{\Delta\text{OD}_S - \Delta\text{OD}_{\text{NADH}}}{\text{OD}_{\text{STD}} - \text{OD}_{\text{BLK}}} \quad (\text{U/L})$$


OD_{STD} and OD_{BLK} are the $\text{OD}_{340\text{nm}}$ values of NADH Standard and Blank at 5 min, respectively. The factor 388 is derived from

$$\text{Factor} = 10 \text{ mM NADH} \times \frac{4 \mu\text{L Vol.}_{\text{NADH}}}{206 \mu\text{L Vol.}_{\text{WR}}} \times \frac{200 \mu\text{L Vol.}_{\text{WR}}}{220 \mu\text{L Vol.}_{\text{Total}}} \times \frac{11 (\text{sample dilution})}{5 \text{ min}} = 388 \mu\text{M/min}$$

If the calculated AST activity is higher than 100 U/L, dilute sample in Assay Buffer and repeat assay. Multiple results by the dilution factor.

MATERIALS REQUIRED, BUT NOT PROVIDED

Pipeting devices and accessories. Clear bottom 96-well plates (e.g. Corning Costar) and plate reader or spectrophotometer and cuvettes for measuring $\text{OD}_{340\text{nm}}$.



Purified AST enzyme (assayed at 37°C)

LITERATURE

- Bergmeyer H.U., Scheibe P. and Wahlefeld A.W. (1978). Optimization of methods for aspartate aminotransferase and alanine aminotransferase. Clin. Chem. 24(1): 58-73.
- Bowers Jr G.N. and McComb R.B. (1984). A unifying reference system for clinical enzymology: aspartate aminotransferase and the International Clinical Enzyme Scale. Clin. Chem. 30(7): 1128-1136.
- Siest G. et al. (1975). Aspartate aminotransferase and alanine aminotransferase activities in plasma: statistical distributions, individual variations, and reference values. Clin Chem. 21(8):1077-87.