

Mouse Pai-1 Total Antigen Assay

Catalog No. IMPAIKT-TOT

Lot No. 612

Intended Use

Murine PAI-1 total assay is intended for the quantitative determination of total plasminogen activator inhibitor type 1 in mouse plasma.

Background

Plasminogen activator inhibitor-1 (PAI-1) is a central regulator of the blood fibrinolytic system [2]. Clinical studies have indicated that increased PAI-1 levels increase the risk for thrombosis, whereas decreased levels may cause recurrent bleeding [3].

Assay Principle

Murine PAI-1 present in plasma reacts with the capture antibody coated and dried on a microtiter plate. Free, latent, and complexed PAI-1 will bind to the plate. Any unbound PAI-1 is washed away and an anti-PAI-1 primary antibody is added. Excess primary antibody is washed away and bound antibody, which is proportional to the total PAI-1 present in the samples, is then reacted with the secondary antibody. Following an additional washing step, TMB is then used for color development at 450nm. The amount of color development is directly proportional to the concentration of total PAI-1 in the sample.

Reagents Provided

Coated plate:

containing mouse anti-mouse PAI-1 antibody, blocked and dried

10X Wash Buffer:

1 bottle of 50ml wash; bring to 1X using DI water

Murine PAI-1 activity standard:

1 vial lyophilized standard

Anti-murine PAI-1 primary antibody:

1 vial lyophilized polyclonal anti-mouse antibody

Anti-rabbit horseradish peroxidase secondary antibody:

1 vial concentrated HRP labeled antibody

TMB substrate solution:

10 ml

Storage and Stability

All kit components must be stored at 4°C. Store unopened plate and any unused microtiter strips in the pouch with desiccant. Reconstituted standards and primary may be stored at -70°C for later use. DO NOT freeze/thaw the standards and primary antibody more than once. All other unused kit components must be stored at 4°C. Kit should be used no later than the expiration date.

Required Reagents and Equipment

1-channel pipettes covering 0-10µL and 200-1000µL
12-channel pipette covering 30-300µL
Paper towels or kimwipes
50mL tubes, 1.5mL centrifuge tubes
1N H₂SO₄
DI water
Magnetic stirrer and stir-bars
Plastic containers with lids
Microtiter plate spectrophotometer operable at 450nm
Microtiter plate shaker with uniform horizontally circular movement up to 300rpm.

Warnings and Precautions

Warning – Avoid skin and eye contact when using TMB substrate solution. It may be irritating to eyes, skin, and respiratory system. Wear safety goggles and gloves.

•DO NOT mix any reagents or components of this kit with any reagents or components of any other kit.

This kit is designed to work properly as provided.

- DO NOT pipette reagents by mouth.
- Always pour substrate out of the bottle into a clean test tube. DO NOT pipette out of the bottle as you could contaminate the substrate.
- Keep plate covered except when adding reagents, washing, or reading.
- DO NOT smoke, drink, or eat in areas where specimens or reagents are being handled.

Preparation of Reagents

TBS buffer: 0.10M TRIS, 0.15M NaCl, pH 7.4

Blocking buffer (BSA): 3% BSA in TBS buffer

Specimen Collection

Collect 9 volumes of blood in 1 volume of a 3.8% trisodium citrate or acidified citrate. Immediately after collection of blood, samples must be centrifuged at 3000Xg for 15 minutes. It is important to ensure a platelet free preparation since platelets can release PAI-1 [4]. The plasma must be stored on ice prior to analysis. The PAI-1 activity samples collected is stable for up to 24 hours or stored at -20 C for up to one month and thawed three times without loss of PAI-1 activity.

Assay Procedure

Perform assay at room temperature. Vigorously shake plate (300rpm) at each step of the assay.

Preparation of Standard:

Reconstitute standard as directed on vial to give a 50ng/ml standard solution.

Dilution table for preparation of mouse PAI standard:

PAI concentration (ng/ml)	Dilutions
50	100µl from standard vial
20	600µl (BB) + 400µl (50ng/ml)
10	500µl (BB) + 500µl (25ng/ml)
5	500µl (BB) + 500µl (10ng/ml)
2.5	500µl (BB) + 500µl (5ng/ml)
1	600µl (BB) + 400µl (2ng/ml)
0.5	500µl (BB) + 500µl (1ng/ml)
0.2	600µl (BB) + 400µl (0.5ng/ml)
0.1	500µl (BB) + 500µl (0.25ng/ml)
0.05	500µl (BB) + 500µl (0.1ng/ml)
0	500µl (BB) Zero point to determine background

Standard and Unknown Addition:

Remove microtiter plate from bag and add 100µl PAI-1 standards (enough for duplicates) and unknowns to wells. Carefully record position of standards and unknowns. Shake plate at 300rpm for 30 minutes. Wash wells three times with 300µl wash buffer. Remove excess wash by gently tapping plate on paper towel or kimwipe.

NOTE: If the unknown is thought to have high PAI-1 levels, dilutions may be made in plasma devoid of PAI-1 or in blocking buffer

Primary Antibody Addition:

Reconstitute primary antibody as directed on vial and agitate gently to completely dissolve contents. Add 100µl to all wells. Shake plate at 300rpm for 30 minutes. Wash wells three times with 300µl wash buffer. Remove excess wash by gently tapping plate on paper towel or kimwipe.

Assay Procedure Continued

Secondary Reagent Addition:

Dilute 1 μ l of conjugated secondary antibody in 10ml of 3% blocking buffer and add 100 μ l to all wells. Shake plate at 300rpm for 30 minutes. Wash wells three times with 300 μ l wash buffer. Remove excess wash by gently tapping plate on paper towel or kimwipe.

Substrate Incubation:

Add 100 μ l TMB substrate to all wells and shake plate for 1-5 minutes. Substrate will change from colorless to different strengths of blue. Quench reaction by adding 50 μ l of 1N H₂SO₄ stop solution to all wells when samples are visually in the same range as the standards. Add stop solution to wells in the same order as substrate upon which color will change from blue to yellow. Mix thoroughly and read final absorbance values at 450nm. For best results read plate immediately.



Measurement:

Set the absorbance at 450nm in a microtiter plate spectrophotometer. Measure the absorbance in all wells at 450nm. Subtract zero point from all standards and unknowns to determine corrected absorbance (A₄₅₀).

Assay Calibration:

Plot A₄₅₀ against the amount of PAI-1 in the standards. Fit a straight line through the points using a linear fit procedure. The PAI-1 activity in the unknowns can be determined from this curve.

Standard Curve Examples



Expected Values

The concentration level of PAI-1 antigen in murine plasma was found to be 1.9 +/- 0.6 ng/ml [1].

Abnormalities in PAI-1 levels have been reported in the following condition:

- Endotoxemia: Endotoxin induces a large increase in PAI-1 levels (80-fold) [1].
- Hyperglycemia, hyperinsulinemia, and insulin resistance: Elevated PAI-1 levels in obese and diabetic mice contribute to these metabolic disorders [5,6].
- Vascular thrombosis: Increased PAI-1 levels may contribute to venous thrombosis [7].
- Myocardial Infarction: Increased PAI-1 levels may contribute to myocardial infarction [7].

Performance Characteristics

Sensitivity = 0.032 ng/ml

(calculated by determining the OD of 20 reps of So and 20 reps of the low standard)

Linearity

The slope = 0.9315

Correlation coefficient = 0.9998

Intra Assay Precision

High 2.6%, Medium 5.2%, Low 7.2%

(calculated by running 20 reps of each concentration in an assay)

Disclaimer

This information is believed to be correct but does not claim to be all-inclusive and shall be used only as a guide. The supplier of this kit shall not be held liable for any damage resulting from handling or from contact with the above product.