

Catalase Activity Assay Kit (Catalog #K773-100; 100 reactions; Store kit at -20°C)

I. Introduction:

Catalase (EC 1.11.1.6) is a ubiquitous antioxidant enzyme that is present in nearly all living organisms. It functions to catalyze the decomposition of hydrogen peroxide (H₂O₂) to water and oxygen. BioVision's Catalase Assay Kit provides a highly sensitive, simple, direct and HTSready assay for measuring Catalase activity in biological samples. In the assay, catalase first reacts with H₂O₂ to produce water and oxygen, the unconverted H₂O₂ reacts with OxiRed™ probe to produce a product, which can be measured at 570 nm (Colorimetric method) or at Ex/Em = 535/587 nm (fluorometric method). Catalase activity is reversely proportional to the signal. The kit detects high pU of catalase activity in samples.

II. Kit Contents:

| Components | K773-100 | Cap Code | Part Number |
|---|----------|----------|-------------|
| Catalase Assay Buffer | 25 ml | NM | K773-100-1 |
| OxiRed™ Probe (lyophilized) | 2 vials | Red | K773-100-2 |
| DMSO (anhydrous) | 0.4 ml | Brown | K773-100-3 |
| HRP (lyophilized) | 1 vial | Green | K773-100-4 |
| H ₂ O ₂ (3%; 0.88M) | 25 µl | Yellow | K773-100-5 |
| Stop Solution | 1 ml | White | K773-100-6 |
| Catalase Positive Control (lyophilized) | 1 vial | Blue | K773-100-7 |

III. Storage and Handling:

Store the kit at -20°C, protect from light. Warm the assay buffer to room temperature before use. Briefly centrifuge vials before opening. Read the entire protocol before performing the assay.

IV. Reagent Reconstitution and General Consideration:

OxiRed™ Probe: Dissolve each vial with 110 µl DMSO (provided, need to warm > 20°C to completely melt) prior to use; one vial is sufficient for 50 assays.

HRP Solution: Dissolve with 220 µl Assay Buffer; it is sufficient for 100 assays.

Positive Control Solution: Dissolve positive control vial into 500 µl Assay Buffer. Aliquot 100 µl per vial, store at -20°C.

Note: Warm the Assay Buffer to room temperature before use. Keep samples, HRP and Catalase on ice during the assay. All these components are stable for 2 week at 4°C or for 1 month at -20°C after reconstitution.

V. Catalase Activity Assay:

1. Sample and Positive Control Preparations:

Homogenize 0.1 gram tissues, or 106 Cells, or 0.2 ml Erythrocytes on ice in 0.2 ml cold Assay Buffer; Centrifuge at 10,000 x g for 15 min at 4°C; Collect the supernatant for assay, keep on ice. Liquid samples can be tested directly. Store samples at -80°C to assay later. Add 2 - 78 µl of samples or 1 - 10 µl Positive Control Solution into each well, and adjust volume to total 78 µl with Assay Buffer. Prepare sample High Control (HC) with the same amount of sample in separate wells then bring total volume to 78 µl with Assay Buffer. Add 10 µl of Stop Solution into the sample HC, mix and incubate at 25°C for 5 min to completely inhibit the

catalase activity in samples as High Control. For unknown samples, we suggest testing several doses of your sample to ensure the readings are within the linear range.

Reducing agents in samples may interfere with the assay. Keep DTT or 2-mercaptoethanol below 5 µM.

2. H₂O₂ Standard Curve:

Dilute 5 µl of 0.88M H₂O₂ into 215 µl dH₂O to generate 20 mM H₂O₂, then take 50 µl of the 20

mM H₂O₂ and dilute into 0.95 ml dH₂O to generate 1 mM H₂O₂. Add 0, 2, 4, 6, 8, 10 µl of 1 mM

H₂O₂ solution into 96-well plate to generate 0, 2, 4, 6, 8, 10 nmol/well H₂O₂ standard. Bring the

final volume to 90 µl with Assay Buffer. Add 10 µl Stop Solution into each well. For the fluorometric assay, dilute the standard H₂O₂ 10-fold for the standard curve (0-1 nmol range).

Note: Diluted H₂O₂ is unstable, prepare fresh dilution each time.

3. Catalase Reaction:

Add 12 µl fresh 1 mM H₂O₂ into each well of both samples and sample HC to start the reaction,

incubate at 25°C for 30 min, and then add 10 µl Stop Solution into each sample vial to stop the

reaction (Note: High Control and standard curve wells already contain Stop Solution).

4. 4. Develop Mix: Mix enough reagents for the number of assays to be performed. For each well

prepare a 50 µl Developer Mix containing:

46 µl Assay Buffer

2 µl OxiRed™ Probe

2 µl HRP solution

Add 50 µl of the Developer Mix to each test samples, controls, and standards. Mix well and incubate at 25°C for 10 min. Measure OD 570 nm in a plate reader.

For low amounts of catalase, you can either increase the incubate time prior to adding the Stop

Solution or use the fluorometric method. For the fluorometric method, decrease the 1 mM H₂O₂

amount to 1.5 µl and OxiRed™ Probe to 0.3 µl in the reaction; compensate the volume with Assay Buffer.

5. Calculation: Signal changes by catalase in sample is $\Delta A = A_{HC} - A_{Sample}$. A_{HC} is the reading of

sample High Control, A_{Sample} is the reading of sample in 30 min. Plot the H₂O₂ Standard Curve.

Apply the ΔA to the H₂O₂ standard curve to get B nmol of H₂O₂ decomposed by catalase in 30

min reaction. Catalase activity can be calculated:


$$\text{Catalase Activity} = \frac{B}{30 \times V} \times \text{Sample Dilution Factor} = \text{nmol/min/ml} = \text{mU/mL}$$

Where: **B** is the decomposed H₂O₂ amount from H₂O₂ Standard Curve (in nmol).

V is the pretreated sample volume added into the reaction well (in ml).

30 is the reaction time 30 min.

Unit definition: One unit of catalase is the amount of catalase decomposes 1.0 μmol of H₂O₂ per min at pH 4.5 at 25 °C.



VI. Related Products:

Colorimetric Glutathione Detection Kit Glutathione Kit (GSH, GSSG and Total)

GST Assay Kit Triglyceride Assay Kit

Phosphatase Assay Kit ADP/ATP Ratio Assay Kit

Phosphate Assay Kit NAD(P)/NAD(P)H Quantification Kit

Pyruvate Assay Kit Lactate Assay Kit/ II

Ammonia Assay Kit Glutamate Assay Kit

Glucose Assay Kit Fatty Acid Assay Kit

Ethanol Assay Kit Uric Acid Assay Kit

Glycogen Assay Kit Phosphate Assay Kits

GENERAL TROUBLESHOOTING GUIDE

| Problems | Cause | Solution |
|---|---|---|
| Assay not working | <ul style="list-style-type: none"> • Use of ice-cold assay buffer • Omission of a step in the protocol • Plate read at incorrect wavelength • Use of a different 96-well plate | <ul style="list-style-type: none"> • Assay buffer must be at room temperature • Refer and follow the data sheet precisely • Check the wavelength in the data sheet and the filter settings of the instrument • Fluorescence: Black plates (clear bottoms) ; Luminescence: White plates ; Colorimeters: Clear plates |
| Samples with erratic readings | <ul style="list-style-type: none"> • Use of an incompatible sample type • Samples prepared in a different buffer • Cell/ tissue samples were not completely homogenized • Samples used after multiple free-thaw cycles • Presence of interfering substance in the sample • Use of old or inappropriately stored samples | <ul style="list-style-type: none"> • Refer data sheet for details about incompatible samples • Use the assay buffer provided in the kit or refer data sheet for instructions • Use Dounce homogenizer (increase the number of strokes); observe for lysis under microscope • Aliquot and freeze samples if needed to use multiple times • Troubleshoot if needed • Use fresh samples or store at correct temperatures until use |
| Lower/ Higher readings in Samples and Standards | <ul style="list-style-type: none"> • Improperly thawed components • Use of expired kit or improperly stored reagents • Allowing the reagents to sit for extended times on ice • Incorrect incubation times or temperatures • Incorrect volumes used | <ul style="list-style-type: none"> • Thaw all components completely and mix gently before use • Always check the expiry date and store the components appropriately • Always thaw and prepare fresh reaction mix before use • Refer datasheet & verify correct incubation times and temperatures • Use calibrated pipettes and aliquot correctly |
| Readings do not follow a linear pattern for Standard curve | <ul style="list-style-type: none"> • Use of partially thawed components • Pipetting errors in the reaction mix • Air bubbles formed in well • Calculation errors • Substituting reagents from older kits/ lots | <ul style="list-style-type: none"> • Thaw and resuspend all components before preparing the reaction mix • Prepare a master reaction mix whenever possible • Pipette gently against the wall of the tubes • Recheck calculations after referring the data sheet • Use fresh components from the same kit |
| Unanticipated results | <ul style="list-style-type: none"> • Measured at incorrect wavelength • Samples contain interfering substances • Use of incompatible sample type • Sample readings above/below the linear range | <ul style="list-style-type: none"> • Check the equipment and the filter setting • Troubleshoot if it interferes with the kit • Refer data sheet to check if sample is compatible with the kit or optimization is needed • Concentrate/ Dilute sample so as to be in the linear range |
| <p>Note# The most probable list of causes is under each problem section. Causes/ Solutions may overlap with other problems.</p> | | |